

## Syringe Strainer 20 µm





| Order No.   | Sterility   | Variant             |
|-------------|-------------|---------------------|
| 43-71020-50 | Non-Sterile | 12 pcs, non-sterile |
| 43-71020-51 | Sterile     | 12 pcs sterile      |

### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

|           |  |
|-----------|--|
| Mesh Size | 20 µm  |
| 大小        | 12 pcs.  |
| 流量控制      | Yes  |
| 外壳材料      | LD-PE (Low Density Polyethylen)                              |
| Mesh Type | woven  |
| 不育症       | Sterile  |
| 运输条件      | Room Temperature   |
| 储存条件      | Room Temperature   |
| 规范性声明     | For research use only. Not for use in diagnostic procedures. |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

