

ÜberStrainer 100 µm



Order No.:43-70100-03

Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

Additional Information

| | |
|-----------------|---|
| Mesh Size | 100 µm |
| 大小 | 10 pcs |
| Fabric material | PET (Polyethylenterephthalat) |
| 外壳材料 | PP (Polypropylene) |
| Mesh Type | woven |
| 不育症 | Sterile |
| 储存条件 | Room Temperature |
| 规范性声明 | For research use only. Not for use in diagnostic procedures., |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

