

## Syringe Strainer 300 µm





Order No.	Sterility	Variant
43-71300-50	Non-Sterile	12 pcs, non-sterile
43-71300-51	Sterile	12 pcs sterile

### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

Mesh Size	300 µm
大小	12 pcs.
流量控制	Yes
外壳材料	LD-PE (Low Density Polyethylen)
Mesh Type	woven
不育症	Sterile
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

