

## Monocytes Spin Medium, 500 ml

**Order No.:**60-00095-12



### Description

Monocytes Spin Medium is a ready-to-use, sterile density medium for the enrichment of monocytes in a high yield from fresh peripheral blood or buffy coat through single step density centrifugation. Depending on the sample a purity of 60-80% is possible. Monocytes Spin is a gentle alternative to bead-mediated monocyte enrichments, especially for downstream applications like cell culture and differentiation.

### Additional Information

大小	500 ml
交付时间	1-2
储存条件	Room temperature, protect from light
规范性声明	For research use only. Not for use in diagnostic procedures.
相位范围	6.8 - 7.4

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.