

Monocytes Spin Medium



Order No.
Size

60-00095-10

100 ml

60-00095-11

250 ml

60-00095-12

500 ml

Description

Monocytes Spin Medium is a ready-to-use, sterile density medium for the enrichment of monocytes in a high yield from fresh peripheral blood or buffy coat through single step density centrifugation. Depending on the sample a purity of 60-80% is possible. Monocytes Spin is a gentle alternative to bead-mediated monocyte enrichments, especially for downstream applications like cell culture and differentiation.

Additional Information

交付时间	1-2
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
储存条件	Room temperature, protect from light
规范性声明	For research use only. Not for use in diagnostic procedures.
相位范围	6.8 - 7.4

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.