

High Density Spin Medium 500 ml (1.1 g/ml)

Order No.:60-01103-12



Description

Density Diluent Medium is a ready to use, sterile medium for the preparation of different density working solutions for the enrichment of cells, bacteria and other particles. The working solutions can be prepared by the dilution of High Density Spin Medium with Density Diluent Medium.

Additional Information

大小	500 ml
交付时间	1-2
运输条件	Room Temperature
储存条件	Room temperature, protect from light
规范性声明	For research use only. Not for use in diagnostic procedures.
相位范围	6.8 – 7.4

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.