

pluriStrainer Mini 200 µm (Cell Strainer)



	Order No.	Size	Sterility	Packaging
<input type="checkbox"/>	43-10200-40	25 pcs.	Sterile	25 pcs. (25 pcs. / Bag) - The tube for flow cytometry is not included.
<input type="checkbox"/>	43-10200-60	100 pcs.	Sterile	100 pcs. (25 pcs. / Bag) - The tube for flow cytometry is not included.
<input type="checkbox"/>	43-10200-50	50 pcs.	Sterile	50 pcs. (25 pcs. / Bag) - The tube for flow cytometry is not included.
<input type="checkbox"/>	43-10200-70	500 pcs.	Sterile	500 pcs. (25 pcs. / Bag) - The tube for flow cytometry is not included.

Description

The pluriStrainer® Mini is a sieving device (or cell strainer) for a wide range of laboratory applications where filtration and purification of liquids is required. It has been optimized for small sample volumes. Its unique design allows fitting on a wide variation of tubes, e.g. 1.5 ml /2.0 ml reaction tubes, 15 ml conical centrifuge tubes, FACS™ tubes, Cryo vials, 24 well plate and 48 well plate. In combination with a standard flow cytometry tube the pluriStrainer Mini 5µm can be used as replacement for a strainer cap.

Application

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marcol-sm-12, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry(FACS™)

Features

The small size laboratory strainer fits into multiple tube sizes



The picture shows, from the left to right: 15 ml conical centrifuge tubes, blood collecting tube, flow cytometry (FACS™) tube, 3 ml tube, 2 ml cryo tube, 5 ml conical tube screw cap, 5 ml reaction tube snap cap, 2 ml reaction tube, 1.5 reaction tube, 48 well plate, 24 well plate

Additional Information

交付时间	1-2
Mesh Size	200 µm
包装	50 pcs. (25 pcs. / Bag) - The tube for flow cytometry is not included.
颜色	light transparent brown
流量控制	No
Mesh Fixation	Injected in housing for strong hold
Fabric material	PET (Polyethylenterephthalat)
外壳材料	PP (Polypropylene)
Mesh Type	woven
符合	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
处理	Gripping surface enables easy and aseptic handling.
不育症	Sterile
相当于	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
稳定性	RLT® Buffer, TRIzol®, isopropanol, organic solvents
离心	Yes
组织解离	Yes
可消毒	Yes
运输条件	Room temperature
储存条件	Room temperature
规范性声明	For research use only. Not for use in diagnostic procedures

法律信息

Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.