

## Re-Strainer (In-Line-Strainer) 300 µm - Sterile

**Order No.:**43-75300-40



### Description

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

### Application

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

### Additional Information

交付时间	1-2
Mesh Size	300 µm
包装	4 devices per bag
大小	20 pcs.

颜色	red
流量控制	yes
Mesh Fixation	Injected in housing for strong hold
Fabric material	PET (Polyethylenterephthalat)
外壳材料	LD-PE (Low Density Polyethylen)
不育症	Sterile
稳定性	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only
法律信息	<b>The product is for research and development only, not for diagnostic or therapeutic use</b>

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.