

pluriSpin® Human CD4⁺ T Cell Enrichment (10 ml)

Order No.:19-03001-32



Description

pluriSpin® Human CD4⁺ T Cell Enrichment is designed to isolate pure, viable and untouched CD4⁺ T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, cytotoxic T cells, B cells, NK cells and granulocytes, while the negatively isolated human CD4⁺ T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD4⁺ T cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

大小	10 ml
品种	Human
细胞类型	Helper T cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

