

## pluriSpin® Human CD8<sup>+</sup> Cell Depletion (10 ml)

**Order No.:**19-04002-32



### Description

pluriSpin® Human CD8<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched cytotoxic t cell-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Cytotoxic T cells, while the isolated human cytotoxic t cell-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified cytotoxic t cell-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

大小	10 ml
品种	Human
细胞类型	Cytotoxic T cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

