

pluriSpin® Human NK Cell Enrichment (10 ml)

Order No.:19-05001-32



Description

pluriSpin® Human NK Cell Enrichment is designed to isolate pure, viable and untouched NK cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, T cells, B cells and granulocytes, while the negatively isolated human NK cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified NK cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

大小	10 ml
品种	Human
细胞类型	NK cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>80%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

