

Syringe Strainer 1 µm (12 pcs., sterile)

Order No.:43-71001-51



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

| | |
|----------------------|--|
| Mesh Size | 1 µm |
| Variation | 12 pcs sterile |
| Size | 12 pcs. |
| Flow Control | Yes |
| Housing Material | LD-PE (Low Density Polyethylen) |
| Mesh Type | woven |
| Sterility | Sterile |
| Shipping Condition | Room Temperature |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.