

## pluriSpin® Human B Cell Enrichment (2 ml)

**Order No.:**19-06001-31



### Description

pluriSpin® Human B Cell Enrichment is designed to isolate pure, viable and untouched B cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, T cells, NK cells and granulocytes, while the negatively isolated human B cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified B cells are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

|                      |   |
|----------------------|---|
| Size                 | 2 ml  |
| Species              | Human   |
| Cell Types           | B cells   |
| Sample Material      | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity               | >80%  |
| Storage Condition    | 4-8 °C  |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

