

10x PBS Solution (500 ml)

Order No.:60-00010-11



Descrizione

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- RNase/DNase free
- Protease free

Components

- NaCl (136,89 mMol/L)
- Na₂HPO₄ • 2H₂O (10mMol/L)
- KH₂PO₄ (1,98 mMol/L)
- KCl (2,68 mMol/L)

Additional Information

Dimensione	500 ml
Sterilità	Sterile
Gamma di ph	pH 6.8
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	Temperatura ambiente

Dichiarazione di regolamentazione

For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.