

BSA stock solution (50 ml)

Order No.:60-00020-10



Descrizione

10% BSA stock solution in 1x PBS can used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

Additional Information

Dimensione	50 ml
Tempi di consegna	1-2
Area di interesse	Biologia cellulare
Sterilità	Sterile
Gamma di ph	pH 7.4
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	4-8 °C
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.