

EDTA Solution (10 ml)

Order No.:60-00030-10



Descrizione

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as an chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

Additional Information

Dimensione	10 ml
Area di interesse	Biologia cellulare
Sterilità	Sterile
Gamma di ph	pH 7.4
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.