

## pluriPlix



**Order No.:**50-01010-80

### Descrizione

Universal rotation mixer without own drive system for usage on top of a magnetic stirrer (not included in delivery). Particularly designed for the gentle mixing of cells. Allows for 2 mixing options: over-head and tilting-rotation (plus: adjustable rotation angle).

### Additional Information

Dimensione	1 device
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.