

1x PBS Solution (500 ml)



Order No.:60-00000-10

Descrizione

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

Components

- NaCl (136,89 mMol/L)
- Na₂ HPO₄ • 2H₂ O (10mMol/L)
- KH₂ PO₄ (1,98 mMol/L)
- KCl (2,68 mMol/L)

Additional Information

Dimensione	500 ml
Tempi di consegna	1-2
Area di interesse	Biologia cellulare
Sterilità	Sterile
Gamma di ph	pH 7.4
Condizioni di spedizione	Temperatura ambiente

Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.