

S-pluriBead® Maxi Reagent Kit

Order No.:70-50010-12



Descrizione

Maxi S-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer, detachment activation buffer, 10x 1,8ml detachment buffer concentrate, 10 S-pluriStrainer, 10 Funnel, 10 connector rings

Additional Information

Dimensione	10
Tempi di consegna	1-2
Specie	Umano
Materiale campione	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purezza	>90%
Condizione di conservazione	Room Temperature; Buffer 4-8 °C; Buffer D -20°C
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.