

## universal S-pluriBead® anti-ms

**Order No.:**31-GTaMS-10



### Descrizione

pluriBead universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

### Additional Information

Dimensione	2 ml
Specie	Mouse (Murine)
Tipi di celle	Qualsiasi
Materiale campione	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purezza	>90%
Condizione di conservazione	4-8 °C
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.