


## 10x PBS Solution





Order No.	Size	Sterility
60-00010-10	100 ml	Sterile
60-00010-11	500 ml	Sterile

### Descrizione

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without**:

- Calcium
- Magnesium
- Phenol Red
- Protease free

### Components

- NaCl (1368,9 mMol/L)
- Na<sub>2</sub> HPO<sub>4</sub> • 2H<sub>2</sub> O (100mMol/L)
- KH<sub>2</sub> PO<sub>4</sub> (19,8 mMol/L)
- KCl (26,8 mMol/L)

### Additional Information

Tempi di consegna	1-2
Area di interesse	Biologia cellulare
Sterilità	Sterile
Gamma di ph	pH 6.8
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	Temperatura ambiente

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.