

## pluriSpin® Human CD19<sup>+</sup> Cell Depletion


**Order No.**
**Size**

19-06002-31

2 ml

19-06002-32

10 ml

### Descrizione

pluriSpin® Human CD19<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched b cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels B cells, while the isolated human b cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified b cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Specie	Umano
Tipi di celle	Celle B
Materiale campione	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purezza	>90%
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.