

## pluriSpin® Human CD45<sup>+</sup> Cell Depletion


**Order No.**
**Size**

19-09002-31

2 ml

19-09002-32

10 ml

### Descrizione

pluriSpin® Human CD45<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched CD45<sup>+</sup>-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45<sup>+</sup> cells, while the isolated population free of CD45<sup>+</sup> are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45<sup>+</sup> cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Specie	Umano
Tipi di celle	CD45- Cella
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.