

pluriSpin® Human CD8⁺ T Cell Enrichment


Order No.
Size

19-04001-31

2 ml

19-04001-32

10 ml

Descrizione

pluriSpin® Human CD8⁺ T Cell Enrichment is designed to isolate pure, viable and untouched cytotoxic T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, helper T cells, B cells, NK cells and granulocytes, while the negatively isolated human cytotoxic T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified cytotoxic T cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|-----------------------------------|---|
| Specie | Umano |
| Tipi di celle | Cellule T citotossiche |
| Materiale campione | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purezza | >90% |
| Condizione di conservazione | Temperatura ambiente |
| Dichiarazione di regolamentazione | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.