

pluriSpin® Human Monocyte Depletion


Order No.
Size

19-01002-31

2 ml

19-01002-32

10 ml

Descrizione

pluriSpin® Human Monocyte Depletion is designed to isolate pure, viable and untouched monocytes-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels monocytes, while the isolated human monocytes-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Specie	Umano
Tipi di celle	Monociti
Materiale campione	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purezza	>90%
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.