

pluriSpin® Human Monocyte Enrichment


Order No.
Size

19-01001-31

2 ml

19-01001-32

10 ml

Descrizione

pluriSpin® Human Monocyte Enrichment is designed to isolate pure, viable and untouched monocytes from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, B cells, NK cells and granulocytes, while the negatively isolated human monocytes are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|-----------------------------------|---|
| Specie | Umano |
| Tipi di celle | Monociti |
| Materiale campione | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purezza | >90% |
| Condizione di conservazione | Temperatura ambiente |
| Dichiarazione di regolamentazione | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.