

## pluriSpin® Human T Cell Enrichment





| Order No.   | Size  |
|-------------|-------|
| 19-02001-31 | 2 ml  |
| 19-02001-32 | 10 ml |

### Descrizione

pluriSpin® Human T Cell Enrichment is designed to isolate pure, viable and untouched T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, B cells, NK cells and granulocytes, while the negatively isolated human T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified T cells are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

|                                   |   |
|-----------------------------------|---|
| Specie                            | Umano   |
| Tipi di celle                     | Celle a T   |
| Materiale campione                | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purezza                           | >90%  |
| Condizione di conservazione       | Temperatura ambiente  |
| Dichiarazione di regolamentazione | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.