

Centrifuge Tube 50 ml (sterile, 300 pcs.)



Order No.:05-00001-01

Descrizione

50 ml tubes are suitable for a wide range of laboratory applications such as sample preparation or sample storage. The tubes are ideal for the use with our **pluriStrainers** and guarantee a perfect fit with the Connector Ring. Moreover the tubes can be used as container for **density gradient centrifugation** or the enrichment of primary cells with **pluriBead**.

Screw tube, 50 ml, 114 x 28 mm, pointed bottom, PP, with printed writing field and scaling, with mounted red cap, sterile and pyrogen-free / endotoxin-free, 25 pcs./bag

Applicazione

- Mixing
- Incubation
- Centrifugation
- Storage
- Collection vessel for filtration and purification of liquids

Additional Information

Dimensione	300 pcs.
Pre-compilato	Unfilled
Tempi di consegna	1-2
Variazione	20 pcs, sterile
Imballaggio	25 pcs. / Bag
Condizione di conservazione	Temperatura ambiente
Condizioni di spedizione	Temperatura ambiente

Dichiarazione di regolamentazione

For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.