

Bundle pluriStrainer® Mini 5 µm and Flow Cytometry Tubes



	Order No.	Size	Sterility	Packaging
<input type="checkbox"/>	43-10005-46	25 pcs.	Sterile	25 pcs. (25 pcs. / Bag)
<input type="checkbox"/>	43-10005-56	50 pcs.	Sterile	50 pcs. (25 pcs. / Bag)
<input type="checkbox"/>	43-10005-66	100 pcs.	Sterile	100 pcs. (25 pcs. / Bag)
<input type="checkbox"/>	43-10005-76	500 pcs.	Sterile	500 pcs. (25 pcs. / Bag)

Descrizione

Bundle includes sterile pluriStrainers Mini with 5 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 5 µm polyethyleneterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 5 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

Applicazione

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marrow, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry (FACS™)

Additional Information

Dimensione della maglia	5 µm
Imballaggio	50 pcs. (25 pcs. / Bag)
Colore	giallo
Controllo del flusso	No
Mesh Fixation	Injected in housing for strong hold
Materiale del tessuto	PET (polietilene tereftalato)
Materiale dell'alloggiamento	PP (Polipropilene)

Si inserisce in	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
Manipolazione	Gripping surface enables easy and aseptic handling.
Sterilità	Sterile
Paragonabile con	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
Stabilità	RLT® Buffer, TRIzol®, isopropanol, organic solvents
Centrifugabile	Sì
Dissociazione dei tessuti	Sì
Disinfettabile	Sì
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	Temperatura ambiente
Informazioni legali	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.