

## ÜberStrainer Set 1

**Order No.:**43-70000-97



### Descrizione

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 20 µm, 60 µm, 100 µm, 300 µm, 500 µm.

### Additional Information

|                                   |  |
|-----------------------------------|--|
| Dimensione della maglia           | Multiple   |
| Dimensione                        | 10 pcs   |
| Materiale del tessuto             | PET (polietilentereftalato)                                  |
| Materiale dell'alloggiamento      | PP (Polipropilene)   |
| Mesh Type                         | woven  |
| Sterilità                         | Sterile  |
| Condizioni di spedizione          | Temperatura ambiente   |
| Condizione di conservazione       | Temperatura ambiente   |
| Dichiarazione di regolamentazione | For research use only. Not for use in diagnostic procedures. |

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.