

## ÜberStrainer Set 3

**Order No.:**43-70000-99



### Descrizione

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 100 µm, 200 µm, 300 µm, 400 µm, 500 µm.

### Additional Information

Dimensione della maglia	Multiple
Dimensione	10 pcs
Materiale del tessuto	PET (polietilentereftalato)
Materiale dell'alloggiamento	PP (Polipropilene)
Mesh Type	woven
Sterilità	Sterile
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures.

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.