

## Syringe Strainer 60 µm (12 pcs., non-sterile)

**Order No.:**43-71060-50



### Descrizione

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

Dimensione della maglia	60 µm
Variazione	12 pcs, non-sterile
Dimensione	12 pcs.
Controllo del flusso	Sì
Materiale dell'alloggiamento	LD-PE (Polietilene a bassa densità)
Mesh Type	woven
Sterilità	Non-Sterile
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.