

Syringe Strainer 30 µm



Order No.	Sterility	Variant
43-71030-50	Non-Sterile	12 pcs, non-sterile
43-71030-51	Sterile	12 pcs sterile

Descrizione

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

Dimensione della maglia	30 µm
Dimensione	12 pcs.
Controllo del flusso	Sì
Materiale dell'alloggiamento	LD-PE (Polietilene a bassa densità)
Mesh Type	woven
Sterilità	Sterile
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

