

Re-Strainer (In-Line-Strainer) 1 µm - Sterile

Order No.:43-75001-40



Descrizione

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

Applicazione

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

Additional Information

Tempi di consegna	1-2
Dimensione della maglia	1 µm
Variazione	20 pcs, sterile
Imballaggio	4 devices per bag

Dimensione	20 pcs.
Colore	rosso
Controllo del flusso	sì
Mesh Fixation	Injected in housing for strong hold
Materiale del tessuto	PET (polietilentereftalato)
Materiale dell'alloggiamento	LD-PE (Polietilene a bassa densità)
Sterilità	Sterile
Stabilità	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
Condizioni di spedizione	Room Temperature
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only
Informazioni legali	The product is for research and development only, not for diagnostic or therapeutic use

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.