

## Re-Strainer (In-Line-Strainer) 20 µm - Sterile

**Order No.:**43-75020-40



### Descrizione

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

### Applicazione

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

### Additional Information

Tempi di consegna	1-2
Dimensione della maglia	20 µm
Imballaggio	4 devices per bag
Dimensione	20 pcs.

Colore	verde
Controllo del flusso	sì
Mesh Fixation	Injected in housing for strong hold
Materiale del tessuto	PET (polietilentereftalato)
Materiale dell'alloggiamento	LD-PE (Polietilene a bassa densità)
Sterilità	Sterile
Stabilità	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
Condizioni di spedizione	Room Temperature
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only
Informazioni legali	<b>The product is for research and development only, not for diagnostic or therapeutic use</b>

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.