

## pluriSpin® Human CD8<sup>+</sup> T Cell Enrichment (10 ml)

**Order No.:**19-04001-32



### Descrizione

pluriSpin® Human CD8<sup>+</sup> T Cell Enrichment is designed to isolate pure, viable and untouched cytotoxic T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, helper T cells, B cells, NK cells and granulocytes, while the negatively isolated human cytotoxic T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified cytotoxic T cells are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Dimensione	10 ml
Specie	Umano
Tipi di celle	Cellule T citotossiche
Materiale campione	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purezza	>90%
Condizione di conservazione	4-8 °C
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

