

pluriSpin® Human B Cell Enrichment (2 ml)

Order No.:19-06001-31



Descrizione

pluriSpin® Human B Cell Enrichment is designed to isolate pure, viable and untouched B cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, T cells, NK cells and granulocytes, while the negatively isolated human B cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified B cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Dimensione	2 ml
Specie	Umano
Tipi di celle	Celle B
Materiale campione	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purezza	>80%
Condizione di conservazione	4-8 °C
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

