

## 10x PBS Solution (500 ml)

**Order No.:**60-00010-11



### 概要

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- RNase/DNase free
- Protease free

### Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub>HPO<sub>4</sub> • 2H<sub>2</sub>O (10mMol/L)
- KH<sub>2</sub>PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Additional Information

大きさ	500 ml
不妊症	無菌状態
フォーンレンジ	pH 6.8
うんぱんじょうたい	室温
保管状態	室温

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.