

## BSA stock solution (50 ml)

**Order No.:**60-00020-10



### 概要

10% BSA stock solution in 1x PBS can be used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

## Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

### Additional Information

大きさ	50 ml
納期	1-2
興味領域	細胞生物学
不妊症	無菌状態
フォーンレンジ	pH 7.4
うんぱんじょうたい	室温
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.