

## EDTA Solution (10 ml)

**Order No.:**60-00030-10



### 概要

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as a chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

### Additional Information

大きさ	10 ml
興味領域	細胞生物学
不妊症	無菌状態
フォーレンジ	pH 7.4
うんぱんじょうたい	室温
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.