

## 1x PBS Solution (500 ml)



**Order No.:**60-00000-10

### 概要

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

### Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub>HPO<sub>4</sub> • 2H<sub>2</sub>O (10mMol/L)
- KH<sub>2</sub>PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Additional Information

大きさ	500 ml
納期	1-2
興味領域	細胞生物学
不妊症	無菌状態
フオーンレンジ	pH 7.4
うんぱんじょうたい	室温

保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.