

M-pluriBead® Mini Reagent Kit

Order No.:70-50010-21



概要

Mini M-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer, detachment activation buffer, 5x 1,8ml detachment buffer concentrate, 5 M-pluriStrainer, 5 Funnel, 5 connector rings

Additional Information

大きさ	5
納期	1-2
しゅるい	人間
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	Room Temperature; Buffer 4-8 °C; Buffer D -20°C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.