

## CD11c M-pluriBead® anti-ms

**Order No.:**29-011c0-20



### 概要

CD11c pluriBead® suspension for the positive cell isolation of CD11c+ target cells from any murine sample material. (2ml)

### Additional Information

大きさ	2 ml
しゅるい	マウス
細胞型	活性化T細胞, 樹状細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.