

CD19 M-pluriBead® anti-ms

Order No.:29-01900-20



概要

CD19 pluriBead® suspension for the positive cell isolation of CD19+ target cells from any murine sample material. (2ml)

Additional Information

| | |
|--------|---|
| 大きさ | 2 ml |
| しゅるい | マウス |
| 細胞型 | B細胞 |
| サンプル材料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 純潔 | >90% |
| 保管状態 | 4-8 °C |
| きせいぶん | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.