

CD326 M-pluriBead® anti-ms

Order No.:29-32600-20



概要

CD326 pluriBead® suspension for the positive cell isolation of CD326+ target cells from any murine sample material. (2ml)

Additional Information

大きさ	2 ml
しゅるい	マウス
細胞型	循環腫瘍細胞, CTC (循環腫瘍細胞, 上皮細胞)
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.