

## Ly6G M-pluriBead® anti-ms

**Order No.:**29-Ly6G0-20



### 概要

Ly6G pluriBead® suspension for the positive cell isolation of Ly6G+ target cells from any murine sample material. (2ml)

### Additional Information

大きさ	2 ml
しゅるい	マウス
細胞型	顆粒球
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.