

CD44 M-pluriBead® anti-hu

Order No.:19-04400-20



概要

CD44 pluriBead® suspension for the positive cell isolation of CD44+ target cells from any human sample material. (2ml)

Additional Information

大きさ	2 ml
しゅるい	人間
細胞型	白血球, ほとんどの組織細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.