

CD271 M-pluriBead® anti-hu

Order No.:19-27100-20



概要

CD271 pluriBead® suspension for the positive cell isolation of CD271+ target cells from any human sample material. (2ml)

Additional Information

大きさ	2 ml
しゅるい	人間
細胞型	中胚葉細胞, シュワン細胞, 幹・間質細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.