

CD8 S-pluriBead® anti-ms

Order No.:29-00800-10



概要

CD8 pluriBead® suspension for the positive cell isolation of CD8+ target cells from any murine sample material. (2ml)

Additional Information

大きさ	2 ml
しゅるい	マウス
細胞型	N K細胞 (サブセット), 細胞障害性T細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.