

## CD43 M-pluriBead® anti-hu

**Order No.:**19-04300-20



### 概要

CD43 pluriBead® suspension for the positive cell isolation of CD43+ target cells from any human sample material. (2ml)

### Additional Information

大きさ	2 ml
しゅるい	人間
細胞型	単球, T細胞, 顆粒球
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
うんぱんじょうたい	室温
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.